


MrT 3300



PARAmed
MEDICAL SYSTEMS



Best performance

MrJ3300 redefines the standard in dedicated MRI by raising image quality to superior levels.

- **0.32 Tesla** magnet
- **Multi-array** dedicated coils
- State-of-the-art pulse sequences for MSK imaging
- Digital platform with management of **up to 16 channels**

High productivity

Extremely easy patient positioning and immediate set-up ensure a **high workflow**.
Rapid pulse sequences result in shorter exam times.



Unparalleled comfort

MrJ3300 provides outstanding patient comfort, unequalled in its category.

- **Largest magnet opening**
- **Illuminated gantry** provides a relaxing environment
- **Versatile patient table** - scan right or left limb with same patient orientation
- Ergonomically designed coils and positioning accessories

Claustrophobic reactions are reduced to the minimum because of the patient friendly environment provided.

The next step in musculoskeletal MRI

Featuring **0.32 Tesla magnet**, anatomy specific multi-array coils, and an extra-spacious magnet opening, **MrJ3300** is the **new leader** in Dedicated Musculoskeletal MRI.

- the most **advanced**
- the most **comfortable**
- the most **versatile**
- the most **compact**






**C-Spine
Fast Rise T2**




**L-Spine
Fast Rise T2**




**Hip flex
SE T1**

**Shoulder
SE T1**



**Knee
GFE STIR**

**Hand/Wrist
GFE T1**



Musculoskeletal imaging

- Complete **Orthopedic** applications, including Shoulder and Hip.
- Comprehensive imaging of **Cervical Spine** and **Lumbar Spine**.
- Advanced pulse sequences (3D GBASS, 3D STSS, FWS)

Optimized workflow

- **Flash Localizer System** with real time display
- **LED** centering device
- Full DICOM connectivity (print, store, worklist-MPPS)




Reduced installation space

MrJ3300 is the most compact system in its category. It can be installed in a very small area.

Reduced cost of operation

Total Cost of Ownership is significantly low due to ease of installation and minimal power consumption.





PARAMed
MEDICAL SYSTEMS

Paramed srl

16152 Genova - Italy
Corso F.M. Perrone, 73r
Tel +39 010 7404530 • Fax +39 010 7404530
www.paramedmedicalsistemas.com
info@paramedmedicalsistemas.com

Paramed Medical Systems Inc.

6204 W. Oakton Street, Morton Grove, IL 60053 - USA
Toll Free: +1 866 840-7565 • Fax +1 847 470-0580
www.paramedmedicalsistemas.com
info@paramed-usa.com

